# Digital Water Baths

SWB6D SWB15D SWB24D


Instructions for use

Mode d'emploi

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**SWB6D Digital Water Bath**

**SWB24D Digital Water Bath**

# DIGITAL WATER BATHS

SWB6D, SWB15D and SWB24D

## Instructions for use

### Introduction

Thank you for purchasing this piece of Stuart equipment. To get the best performance from the equipment and for your own safety, please read these instructions carefully before use. Before discarding the packaging check that all parts are present and correct.

This equipment is designed to operate under the following conditions:

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Over-voltage category II IEC60364-4-443
- ❖ Pollution degree 2

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Stuart, the protection provided may be impaired.

### Electrical Installation

**THIS EQUIPMENT MUST BE EARTHED**

**Before connection please ensure that the line supply corresponds to that shown on the rating plate.**

| Model  | Supply requirements | Power |
|--------|---------------------|-------|
| SWB6D  | 230V, 50Hz          | 350W  |
| SWB15D | 230V, 50Hz          | 750W  |
| SWB24D | 230V, 50Hz          | 1400W |

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other. Should neither lead be suitable, take the lead with the U.K. plug and

replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the re-wireable plug in accordance with its instructions.

**THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN**

**NOTE:** Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated

**CAUTION – Fuse fitted in both live and neutral lines.**

#### Fuses

230V: 2 x 10A fuses. For all other voltages always replace with an equivalent fuse type and rating as those being removed, alternatively see rating plate for fuse type and rating.

The wires in the mains cable are coloured as follows:

|              |   |         |
|--------------|---|---------|
| BROWN        | - | LIVE    |
| BLUE         | - | NEUTRAL |
| GREEN/YELLOW | - | EARTH   |

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected **BEFORE** connection to the mains supply. When wishing to isolate the bath from the mains supply, disconnect the cable from the IEC socket.

Should the mains lead need replacement, a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC320 plug should be used.

**IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

# Safety Advice Before Use

**USE CAUTION WHEN REMOVING VESSELS FROM THE BATH AS PARTS MAY BE HOT.**

**BE AWARE THAT STEAM MAY BE RELEASED WHEN REMOVING THE LID**


- ❖ The unit should be carried with both hands, with fingers under each side.
- ❖ Never move or carry the unit when filled with water, while it is in use or connected to the mains supply.
- ❖ Use only water and no other liquid in the bath.
- ❖ Do not lean over the bath when removing the lid as steam may be released. Open the lid slowly to allow steam to disperse safely.
- ❖ Allow water to cool before draining.

## General Description


The SWB series digital water baths are ideal for heating vessels and maintaining samples at a constant temperature. They are available in 6L, 15L and 24L options to accommodate different vessels of different sizes. They have a one-piece stainless steel construction for ease of maintenance and cleaning. In addition to this, a drain tap aids removal of the bath water. For safety, all baths have a water level sensor to prevent the bath running dry. Temperature is maintained by a PRT temperature sensor and there is an over-temperature cut out for added safety. To avoid evaporation and improve temperature control, clear polycarbonate covers are included with the baths.

## Controls


### Control panel

### Bath inside view

### Drain tap (on rear)

## Operation

### Setting up

Position the apparatus on a firm level surface and insert the mains lead into the socket on the rear of the unit. Ensure the isolation switch is easily accessible during use.

Attach a suitable length of hose to the drain tap at the rear of the bath to allow for draining off the water after use. Ensure the drain tap is accessible to prevent spillage when the bath is full.

**Note:** To ensure correct operation, prior to use you should clean the level sensor with a mild solvent, such as Isopropyl alcohol.

Connect to the mains electricity supply.

Fill with water to a level above the low water level sensor inside the bath. It is advisable to use distilled or deionised water with neutral pH to avoid scale formation. Ultra pure deionised or RO water with a resistivity >1Mohm should not be used, as this can leach the iron from the stainless steel and lead to corrosion.

Periodically top up with water to maintain the desired level.

### Temperature setting

When connected to the mains electricity supply the LED display will illuminate and show the word OFF.

Switch the unit on by pressing the "Power On/Off" knob. The display will momentarily show the set temperature followed by the actual water temperature in degrees Celsius.

Press and hold the "Press to set" button. The display will show the set temperature. The set temperature may be adjusted by turning the "Temperature adjustment" knob whilst continuing to hold in the "Press to set" button. When the desired setting is displayed release the "Press to set" button. The display will revert to show the actual water temperature and the bath will begin to heat.

The set temperature can be observed at any time by pressing the "Press to set" button.

Turning the adjustment knob without pressing the "Press to set" button will have no effect on the temperature setting.

### Hot indicator

When the bath temperature rises above 50°C, the hot indicator will begin to flash. This will continue even when the power is turned off at the "Power On/Off" knob until the water temperature falls below 50°C. When the power is off, the LED display will also glow with the word "Hot" until the water temperature is below 50°C.

**CAUTION - Use protective gloves to remove hot vessels.**

### Low water level sensor

The low water level sensor is designed to provide a warning when the water level has become too low in order to protect the bath from overheating. When the water level drops below the sensor, the LED display will flash "H2O" and "ERR" to indicate that the water level is too low and the heater will switch off. The heater may be re-started by filling the bath above the level of the sensor.

### Drain tap

The drain tap allows for easy removal of the bath water for maintenance and cleaning. Ensure the bath is turned off and the water is cool before draining.

Direct the hose from the drain trap to a sink or waste collection vessel. Open the tap valve and allow the water to flow out. Close the tap before re-filling.

## Maintenance & Servicing

**WARNING:** Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

### Maintenance

It is recommended that the bath is drained, cleaned and re-filled on a weekly basis to avoid build up of salts and contamination. For baths that are used for biological applications and heated to physiological temperatures, regular heating to >60°C for 30 minutes can be used to thermally disinfect the bath before cleaning. Do not use bleach (sodium hypochlorite) or other chlorine based solutions as a disinfectant.

### Cleaning

In general baths should be cleaned with a mild household or laboratory detergent using a sponge or soft cloth. Scouring powders, steel wool or other abrasive pads should not be used. Scale build up can be removed using a mild household descaler suitable for stainless steel; follow the manufacturer's instructions for correct use. Rinse thoroughly after cleaning and dry the bath.

For cleaning of the low water level sensor, Freon and alcohol based solvents are recommended. DO NOT use chlorinated solvents such as trichloroethane as these are likely to attack the sensor material.

Periodically clean the outside of the unit using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

### Servicing

**Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.**

For a comprehensive list of parts required by service engineers conducting internal repairs, please contact the Service Department, quoting both the model number and serial number:

E-mail: [service@bibby-scientific.com](mailto:service@bibby-scientific.com)

Tel: +44 (0)1785 810475

Fax: +44 (0)1785 810471

**Only spare parts supplied or specified by Stuart or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.** If in any doubt, please contact the Service Department of Bibby Scientific Ltd. or the point of sale.

For any other technical enquiries, please contact the Technical Support Department:  
E-mail: [stuarthelp@bibby-scientific.com](mailto:stuarthelp@bibby-scientific.com)  
Tel: +44 (0)1785 810433.

## Spares and Accessories

| Description                  | Catalogue Number |
|------------------------------|------------------|
| Polycarbonate lid 6L bath    | SWB6D/1          |
| Polycarbonate lid 15L bath   | SWB15D/1         |
| Polycarbonate lid 24L bath   | SWB24D/1         |
| Stainless steel lid 6L bath  | SWB1/1           |
| Stainless steel lid 15L bath | SWB2/1           |
| Stainless steel lid 24L bath | SWB3/1           |

## Warranty

Stuart warrants this equipment to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of **three (3) years**. In the event of a justified claim, Stuart will replace any defective component or replace the unit free of charge. This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.

Bibby Scientific Ltd.  
Beacon Road,  
Stone, Staffordshire  
ST15 0SA  
United Kingdom  
Tel: +44 (0) 1785 812121  
Fax: +44 (0) 1785 813748  
E-mail: [equipment@bibby-scientific.com](mailto:equipment@bibby-scientific.com)  
[www.stuart-equipment.com](http://www.stuart-equipment.com)

## Technical Specification

|                         | SWB6D                  | SWB15D                 | SWB24D                 |
|-------------------------|------------------------|------------------------|------------------------|
| Nominal capacity        | 6L                     | 15L                    | 24L                    |
| Temperature range       | Ambient +5°C to 99.9°C | Ambient +5°C to 99.9°C | Ambient +5°C to 99.9°C |
| Temperature stability   | +/-0.5°C               | +/-0.5°C               | +/-0.5°C               |
| Low water level control | Yes                    | Yes                    | Yes                    |
| Drain tap               | Yes                    | Yes                    | Yes                    |
| Dimensions:             |                        |                        |                        |
| w                       | 335 (300)mm            | 335 (300)mm            | 335 (300)mm            |
| (internal) d            | 230 (150)mm            | 408 (325)mm            | 584 (500)mm            |
| h                       | 280 (200)mm            | 280 (200)mm            | 280 (200)mm            |
| Working depth max/min   | 160/80mm               | 160/80mm               | 160/80mm               |
| Weight                  | 6.5kg                  | 9.0kg                  | 12.8kg                 |
| Heater power            | 350W                   | 750W                   | 1400W                  |
| Electrical supply       | 230V, 50Hz             | 230V, 50Hz             | 230V, 50Hz             |
| IP rating               | 31                     | 31                     | 31                     |

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these

standards and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

## Declaration of Conformity

### Water Bath, Model SWB6D, SWB15D, SWB24D

These products comply with the requirements of the EU Directives listed below:

**2004/108/EC**      **EMC Directive.**  
**2006/95/EC**      **Low voltage Directive (LVD)**

Compliance with the requirements of these Directives is claimed by meeting the following standards:

**EN 61326-1:2006 (Electrical Equipment for Measurement, Control and Laboratory use).**  
**EN 61010-1: 2010**  
**(Safety Requirements Electrical Equipment for Measurement, Control and Laboratory use)**  
**EN 61010-2-010: 2003 (Particular Requirements for Laboratory Equipment for Heating of Materials).**

**CE mark affixed '12.**

Signed:  (Mr S. Marriott)

Date: 21/5/12

Authority: Technical Manager  
Bibby Scientific Ltd

Bibby Scientific Ltd - Stone - Staffs - ST15 0SA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748  
DofC\_88\_Rev01

# INSPECTION REPORT

**MODELS** SWB6D, SWB15D and SWB24D

## ELECTRICAL SAFETY

1. Earth continuity

2. Insulation

3. Flash test

## FUNCTIONAL

1. Indicators


2. Temperature control

3. Visual acceptance

**QUALITY CONTROL INSPECTOR**

**Bibby Scientific**

### **UK**

Bibby Scientific Ltd.  
Beacon Road, Stone  
Staffordshire ST15 0SA  
United Kingdom  
Tel: +44 (0)1785 812121  
Fax: +44 (0)1785 810405  
e-mail: sales@bibby-scientific.com  
www.bibby-scientific.com

### **North & South America**

Bibby Scientific US Inc.  
t/a Techne Inc.  
3 Terri Lane, Suite 10  
Burlington, NJ 08016 USA  
Toll Free (in NA): 800-225-9243  
Tel: +1 609 589 2560  
Fax: +1 609-589-2571  
e-mail: labproducts@techneusa.com  
www.techneusa.com

### **Middle East**

Bibby Scientific Middle East Ltd.  
PO Box 27842, Engomi 2433  
Nicosia  
Cyprus  
Tel: +357 22 660 423  
Fax: +357 22 660 424  
e-mail: sales@bibby-scientificme.com

### **France**

Bibby Scientific Limited  
Bâtiment Le Deltaparc Parc Silic PN2  
7 rue du Canal  
BP 55437 VILLEPINTE  
95944 ROISSY Charles de Gaulle  
France  
Tel: +33(0)148 63 78 03  
Fax: +33(0)148 63 78 01  
e-mail: ventes@bibby-scientific.com  
www.bibby-scientific.com

### **Italy**

Bibby Scientific Italia Srl  
Via Alcide de Gasperi 56  
20077 Riozzo di Cerro al Lambro  
Milano Italia  
Tel: +39 (0)2 98230679  
Fax: +39 (0)2 98230211  
e-mail: marketing@bibby-scientific.it  
www.bibby-scientific.it

### **Asia**

Bibby Scientific - Singapore  
Prudential Tower, Level 26  
30 Cecil Street  
Singapore 049712  
Tel: +65 6631 2976  
Fax: +44 (0)1785 810405  
e-mail: info@bibby-scientific.com  
www.bibby-scientific.com