

CLINIPET+ pipette series

CLINIPET+ fixed volume pipettes are available in 33 models from 2 μL up to 1,000 μL . The handle and mechanism have been especially designed for fixed volume pipettes. The pipette features reduced pipetting forces and excellent accuracy and precision combined with highly durable, UV resistant and fully autoclavable body.

- Recalibration system allows for quick and accurate pipettes calibration



- Low pipetting forces prevent RSI
- Highly durable, made from PVDF shafts guarantee highest mechanical and chemical resistance



- Unique tip sealing system



Sliding tip ejector button – perfect for right & left handed people. Its unique construction reduces tip ejection forces



Easy removable tip ejector allows working with narrow tubes

Model	Cat. no.	Colour code	Volume (µL)	A (%)	P (%)	Standard tips, non-filtered
-------	----------	-------------	-------------	-------	-------	-----------------------------

CP2*	2202-CP		2	± 3.0	± 1.5	
CP5	2204-CP		5	± 2.0	± 1.0	
CP7*	2206-CP		7	± 1.5	± 0.7	
CP10	2208-CP		10	± 1.2	± 0.5	
CP15*	2210-CP		15	± 1.0	± 0.5	
CP20	2212-CP		20	± 0.9	± 0.4	
CP25	2214-CP		25	± 0.9	± 0.4	
CP30*	2216-CP		30	± 0.9	± 0.3	
CP40*	2218-CP		40	± 0.9	± 0.3	
CP44.7*	2220-CP		44.7	± 0.9	± 0.3	
CP50	2222-CP	●	50	± 0.9	± 0.3	
CP60*	2224-CP		60	± 0.9	± 0.3	
CP70*	2226-CP		70	± 0.9	± 0.3	
CP75*	2228-CP		75	± 0.9	± 0.3	
CP80*	2230-CP		80	± 0.9	± 0.3	
CP90*	2232-CP		90	± 0.8	± 0.3	
CP100	2234-CP		100	± 0.8	± 0.3	
CP120*	2236-CP		120	± 0.8	± 0.3	
CP150*	2238-CP		150	± 0.7	± 0.3	
CP200	2240-CP		200	± 0.6	± 0.2	
CP200A*	2242-CP		200	± 0.6	± 0.2	
CP220*	2244-CP		220	± 0.6	± 0.2	
CP250	2246-CP		250	± 0.6	± 0.2	
CP300	2248-CP		300	± 0.6	± 0.2	
CP400*	2250-CP		400	± 0.6	± 0.2	
CP450*	2252-CP		450	± 0.6	± 0.2	
CP500	2254-CP	●	500	± 0.6	± 0.2	
CP600*	2256-CP		600	± 0.6	± 0.2	
CP700*	2258-CP		700	± 0.6	± 0.2	
CP750*	2260-CP		750	± 0.6	± 0.2	
CP800*	2262-CP		800	± 0.6	± 0.2	
CP900*	2264-CP		900	± 0.6	± 0.2	
CP1000	2266-CP		1,000	± 0.6	± 0.2	



200 µL



1,000 µL

* Special version upon request. Note: The HTL maximum permissible errors are guaranteed only when HTL pipette tips are used.



5484

Carousel stand



5449

Multiple stand



5480

Plexi 4-position stand



19451

Pipette hanger