Horizontal cylindrical pressure steam sterilizer

(microcomputer control)



Product details:

Characteristics:

- 1. Sterilizing course automatic sterilization controlled by computer, easy to operate.
- 2. With drying function, suitable dressing drying.
- 3. With over-temperature, over-pressure auto-protect device.
- 4.The door opening mechanism cannot be operated until the pressure in the chamber is reduced to 0.027MPa.It can't be start on if the doesn't closed well.
- 5. The safety value will be open automatically when the inner pressure over 0.24 MPa, and the steam be exhausting to the water tank.
- 6. The power will be cut-off automatically, cut-off the water and alarm when the machine lacking water.
- 7. The chamber of the sterilizer is made of stainless steel.
- 8. With printer.

Technical data:

Technical data	WS-150YDD	WS-200YDD	WS-280YDD
Sterilizing chamber	150L	200L	280L
volume	$\Phi440\!\times\!1000\text{mm}$	ф 515×1000mm	ф 600×1000mm
working pressure	0.22 MPa		
working temperature	134℃		
Adjustment of temperature	40℃ - 134℃		
Time of sterilization	0-99 min		
Time of drying	0-99 min		
Heat average	$\leq \pm 2^{\circ}$		
Power	9KW/380V 50Hz	9KW/380V 50Hz	12KW/380V 50Hz
Dimension	$1400 \times 650 \times 1600$	$1400 \times 750 \times 1700$	$1520 \times 910 \times 1900$
Transport dimension	$1560 \times 820 \times 1820$	$1560 \times 910 \times 1880$	$1680 \times 1080 \times 2100$
G.W/N.W	430/340 KG	436/350 KG	570/462 KG